





TISCH R1002 C BRAVO 2 MIT STAUFACH



EN - 1865



TECHNISCHE DATEN / TECHNICAL DATA

LÄNGE	LENGTH	2089 mm
BREITE	WIDTH	632 mm
HÖHE	HEIGHT	451 mm
GEWICHT	WEIGHT	120 kg
BELADUNG MAX.	MAX. LOAD	370 kg
H. max. Fahrzeugboden	H. max. Ambulance floor	650 mm

EIGENSCHAFTEN / FEATURES

- Tragentisch für TG1000 powerBRAVA
- Aus V2A-Edelstahl und Aluminium
- Automatische Aufnahme der powerBRAVA per Knopfdruck - Easy Loading System / 3-Punkt-Tragenarretierung
- Automatisches Aufladen der Akkus durch Andockstation
- Mit Parallelverschiebung von bis zu 122 mm nach links oder bis zu 132 mm nach rechts
- Tragentisch verfügt über ein Staufach für eine Schaufeltrage oder ein Spine-Board
- Funktion Höhenverstellung der Fahrtrage bei stehendem Fahrzeug
- Tragentisch wird am Fahrzeugboden mit 8 Schrauben fixiert (DIN 7991 M10 - 10.1)



- Rail only for TG-1000 POWER BRAVA stretcher
- Movement toward the exterior of the ambulance to easy loading the stretcher
- Manual displacement on 2 guide lines
- Sideways movement of 254 mm of the whole assembly
- The rail is prepares for both scoop or spinal board
- Rail made of aluminium and stainless steel
- 3 solid anchoring points for the stretcher
- Easy fixation of the rail into the ambulance with an aluminium floor guide

